**BR** **Mettler-Toledo Ind. e Com. Ltda.**, Alameda Araguaia, 451 - Alphaville, BR-06455-000 Barueri/SP, Brazil, Phone +55 11 4166 74 00, Fax +55 11 4166 74 01

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## 1. Introduction

METTLER TOLEDO pH sensors type InPro 3250 are low-maintenance electrodes with pressurized liquid electrolyte, silver-ion trap and built-in temperature sensor. The InPro 3250 sensors are suitable for the same applications as the corresponding DPA(S) sensors without temperature sensor, and provide an additional temperature signal for automatic temperature compensation of the pH signal during calibration and operation. InPro 3250 sensors are autoclavable and in-situ sterilizable. We recommend that the sensor may only be in operation with the original accessories from METTLER TOLEDO. The sensor may be operated and serviced only by personnel familiar with the unit, and who have read and understood this instruction manual.

## 2. Safety remarks

The sensor is intended only for the use described in § 1 above. The sensor consists partly of glass and should be handled carefully. If the sensor is to be cleaned or calibrated using acid or alkaline solutions, safety goggles and gloves should be worn. We recommend that the electrodes only be put into operation in conjunction with original products/accessories from METTLER TOLEDO. The electrodes may be operated and serviced only by personnel familiar with the equipment in question and who have read and understood these instructions. Local safety regulations are applicable to the handling of the sensor. Please follow the guidelines enclosed for an installation in the Ex-zone. Please take extra care during installation/fitting of the sensor to ensure that neither the polymer body/shaft of the sensor nor the thread (Pg 13.5) suffer any injury or damage. Please do not touch or hold the electrode directly at the plug. Electrostatic discharge can damage the chip.

## 3. Product description

The markings on each sensor provide the following information:

<b>METTLER TOLEDO</b>	sensor manufacturer
InPro® 3250 SG/XXX/YYY/ZZZ	designation (Solution Ground Option)/shaft length in mm/temperature sensor/ISM (option)
combination pH	type of sensor
pH 0 ... 14 (InPro 3250), pH 0 ... 12 (InPro 3253)	pH measuring range
0 ... 100 °C (140 °C)	temperature range for operation/sterilization
Order No. 52 00X XXX	article number for ordering
⊗ II 1/2G EEx ia IIC T6/T5/T4/T3	Ex-classification
SNCH 00 ATEX 3130 X	designation and number of certificate
IS/I, II, III/1/ABCEFG/T6	FM approved
Follow the operating instructions	note

Additionally, the sensor is supplied with a serial number on the red connector part for identification possibilities of each sensor produced.

## 4. Installation and preparation for use

1. On unpacking, check the sensor for possible mechanical damage of the pH sensitive glass, the diaphragm, the shaft and the VarioPin (VP) connector. Report any damage immediately to your METTLER TOLEDO supplier.
2. Remove the watering cap and briefly rinse the sensor with de-ionized water. Carefully remove the protection droplet of silicone on the diaphragm with the enclosed knife, as shown on the green pictures on the packing.
3. Again rinse the sensor well with de-ionized water. After rinsing, the sensor should only be dabbed dry with a tissue. Do not rub the pH-sensitive glass, since this can lead to electrostatic charging and sluggish response times.
4. **Carefully check the space behind the pH-sensitive glass for the presence of any air bubbles** and remove same by gently swinging the sensor in a vertical plane.
5. Please refer to the instruction manual of the relative housing in order to install the sensor correctly.
6. Connect the sensor to the pH transmitter, using a cable with VarioPin (VP) connector. Please refer to the wiring and terminal diagram accompanying the cable.
7. For electrodes with Pt100 temperature sensors, 3-wire technique is recommended in order to obtain best measurement accuracy. Please refer to the cable and transmitter instruction manuals for the corresponding cable wire connections.
8. For InPro 3250 SG electrodes (with an additional Pt-ring for solution ground) connect the corresponding blue wire in the VP cable to your transmitter in accordance with instructions in the connection diagram of VP cable and/or in the transmitter instruction manual. For ORP measurement (sep. transmitter needed) connect solution ground with ORP input of the transmitter.
9. InPro 3250 SG electrodes with ISM allow "Plug and Measure" as well as advanced diagnostics. Operation of the electrodes in ISM mode requires the use of a special VP8 cable as well as an ISM pH-module for the transmitter M 700 (software). For installation, startup and operation of the ISM system please refer to the relative instructions provided with the cable and the module.

## 5. Operation

### 5.1 Calibration of the sensor and pH transmitter

For InPro 3250 sensors, a 2-point calibration is recommended. After removing the watering cap, the sensor is dipped in succession into two different buffer solutions with given pH reference values and the transmitter is calibrated to these buffer values.

### 5.2 Sterilization of sensor

In the case of sterile processes, such as in fermentation applications, the sensors are sterilized prior to the production cycle. Autoclaving of the whole fermenter/reactor unit including pH electrode is possible with this sensor type. Sterilization can also be effected in-situ by means of steam or superheated process medium in contact with the lower sensor parts.

### 5.3 Temperature sensor

The integrated Pt100/1000 temperature sensor is primarily intended for automatic temperature compensation of the pH signal during calibration and normal operation, and **not** for the accurate indication or control of the process temperature.

## 6. Maintenance

1. Carefully rinse the sensor tip and diaphragm with de-ionized water after each production cycle. Under no circumstances must measuring solution be allowed to dry on these parts!
2. When the sensor is not in operation, store it with sensor tip and diaphragm well submerged into 9823/3M KCl or 9816/Viscolyt electrolyte.
3. If the sensor is stored mounted in its housing, the same rule as in 2) above applies, but the storage electrolyte has to be slightly modified, with 2 parts of buffer solution pH 9.2 to 10 parts of electrolyte, this in order to avoid corrosion of housing parts.
4. If a sensor is stored dry for a few days by mistake, let it soak in the normal storage electrolyte for several hours before use.
5. Occasionally check the connector for possible traces of moisture. If necessary, clean well with de-ionized water or alcohol, and afterwards dry carefully.

## 7. Trouble-shooting (sluggish response, increased resistance etc.)

**At sluggish response time** (e.g. after a long and dry storage period), reactivation of the pH sensitive glass may be necessary. Dip the pH sensitive glass (but no further sensor parts!) into Reconditioning (HF) Solution (article no. 51319153) for 5 min. Afterwards, rinse well with de-ionized water and let sensor tip and diaphragm soak in reference electrolyte for minimum 12 hours.

**Rapidly increased resistance of reference electrode and/or drift** could be signs of protein precipitation. Use Cleaning Solution Pepsin/HCl (article no. 51 340 068) to remove such precipitation. At repeated problems of this type, consider switching to sensors with Friscolyt™ based gel electrolyte.

## 8. Disposal

If no relative official regulations are in force, used or defective sensors and the packaging can be disposed of as ordinary household waste.

## 9. Warranty

In the event of manufacturing faults, 12 months warranty from date of production is granted.

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