English:

- 1) VP connector
- 2) Pg 13.5 thread
- 3) Washer PTFE
- 4) O-ring (Viton®)
- 5) PEEK (FDA) shaft
- 6) Ceramic diaphragm
- 7) ISFET

Deutsch:

- 1) VP Steckverbindung
- 2) Pg 13.5 Gewinde
- 3) Scheibe PTFE
- 4) O-Ring (Viton®)
- 5) PEEK (FDA) Schaft
- 6) Keramik-Diaphragma
- 7) ISFET

Français:

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Mettler-Toledo GmbH, Process Analytics, Industrie Nord, CH-8902 Urdorf, Tel. +41 1 736 22 11, Fax +41 1 736 26 36
 Subject to technical changes. 07/03. © Mettler-Toledo GmbH. Printed in Switzerland

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General instructions for non-glass pH combination electrodes InPro® 3300

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1. Introduction

METTLER TOLEDO pH sensor type InPro® 3300 is a low-maintenance, autoclavable and pressure-resistant electrode with a glass-free sensing head and thickened reference electrolyte. Please read these instructions carefully before putting the electrode into operation, in order to ensure perfect operation. All wetted parts are made of materials in compliance with the FDA directives (Food and Drug Administration). We recommend that the sensor may only be in operation with the original accessories from METTLER TOLEDO. The sensor may be operated and serviced only by personnel familiar with the unit, and who have read and understood this instruction manual.

2. Safety remarks

If an electrode is to be cleaned or calibrated using acid or alkaline solutions, safety goggles and gloves should be worn.

3. Product designation

The markings on the polymer shaft include the following information:

METTLER TOLEDO	electrode/sensor manufacturer
InPro® 3300/120/Pt xxx	designation / shaft length in mm / temperature sensor
Combination pH	type of electrode
pH 0...14	pH measurement range
0...80 °C / (130 °C)	temperature range for operation/(sterilization)
Order No. 52 xxx xxx	article number for ordering

Additionally, the electrode is supplied with a serial number marked on the connector cap to enable individual identification of each electrode produced.

4. Installation and preparation for use

1. On unpacking, check the electrode for mechanical damage. Report any signs of damage immediately to your METTLER TOLEDO supplier.
2. Remove the watering cap and briefly rinse the electrode with de-ionized water. An already used electrode should first be cleaned by lightly scrubbing the sensing element with a soft-bristled toothbrush and a mild soap solution, followed by rinsing with de-ionized water.
3. Please refer to the instruction manual of the respective housing in order to install the electrode. The electrode is of asymmetrical design. A marking on the connector head shows the direction of alignment of the sensing element for proper installation into the housing.
4. Connect the electrode to the pH transmitter, using cable VP6-ST/ISFET. Please refer to the wiring and terminal diagram accompanying the cable.
5. Prior to calibration, soak the electrode in buffer pH 7.00 or pH 4.01 for at least 10 minutes, with the electrode connected to the transmitter.

5. Operation

5.1 Calibration of the electrode and the pH transmitter

A 2-point pH calibration is recommended e.g. using buffers pH 7.00 and pH 4.01. Please refer to the instruction manual of the pH transmitter for further details.

5.2 CIP-cleanability (Cleaning-In-Place)

A standard procedure for maintaining hygienic conditions in food processing is Cleaning-In-Place (CIP). Periodic flushing with hot caustic soda (NaOH) will however shorten the lifetime of the electrode. We therefore recommend to retract the electrode from the process and to clean it with hot steam (sterilization procedure).

6. Maintenance

1. Carefully rinse the sensor tip and diaphragm with de-ionized water after each production cycle. Under no circumstances must measuring solution be allowed to dry on these parts!
2. When the sensor is not in operation, store it with sensor tip and diaphragm well submerged into 9823/3M KCl or 9816/Viscolyt electrolyte.
3. If the sensor is stored mounted in its housing, the same rule as in 2) above applies, but the storage electrolyte has to be slightly modified, with 2 parts of buffer solution pH 9.2 to 10 parts of electrolyte, this in order to avoid corrosion of housing parts.
4. If a sensor is stored dry for a few days by mistake, let it soak in the normal storage electrolyte for several hours before use.
5. Occasionally check the connector for possible traces of moisture. If necessary, clean well with de-ionized water or alcohol, and afterwards dry carefully.

7. Trouble-shooting

For cleaning, lightly scrub the sensing element of the electrode with a soft toothbrush and water to remove any possible residues. Use mild detergent if necessary. Possible protein precipitation and blockage of the diaphragm can be treated with cleaner 51 340 070.

Attention: HF-containing solutions will destroy the pH sensitive element of the electrode.

8. Disposal

If no official regulations are in force, used or defective electrodes and the packaging can be disposed as ordinary household waste.

9. Warranty

In the event of manufacturing faults, a warranty period of 12 months from date of production is granted.

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